

TRACHEOTEC® VARIO CUFF-SUCTION

- Tracheostomy tube with cuff and integrated suction line to remove secretions above the cuff
- Suction line integrated into the inner wall of the tube to avoid any additional irritation of the trachea or the areas around the tracheostoma
- Soft, flexible neck-flange with key tube identification data clearly marked
- ➤ With rotating 15 mm connector (VARIO) to buffer torque and tension caused by tractive forces from ventilation equipment
- With X-ray contrast stripe
- Set includes tube holder, obturator and product ID card



Order information:



TRACHEOTEC® VARIO CUFF-SUCTION SIZE 5-10

When ordering, please state the size after the item number!

REF 18930

STERILE EO	Sterilized using ethylene oxid
8	Do not resterilize
*	Keep away from sunlight
SANGATA PATECATA USE 1	Single patient use
[]i	Consult instructions for use
Rx only	CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner (in the US only).

CE marking identification number of the notified body

Size	O.D. neck flange mm	O.D. Tip mm	I.D. mm	Length mm (Midline)	Length mm	REF
5	7.3	7.3	5.0	55.7	61.1	18930-05
5.5	8.0	8.0	5.5	57.7	63.6	18930-055
6	8.7	8.7	6.0	60.7	97.1	18930-06
6.5	9.3	9.3	6.5	66.6	73.5	18930-065
7	10.0	10.0	7.0	69.6	77.0	18930-07
7.5	10.7	10.7	7.5	73.6	81.5	18930-075
8	11.3	11.3	8.0	79.5	87.9	18930-08
8.5	12.0	12.0	8.5	81.5	90.4	18930-085
9	12.7	12.7	9.0	85.5	94.9	18930-09
10	14.0	14.0	10.0	95.0	105.3	18930-10

O.D. Shield = outer diameter behind the neck flange · O.D. Tip = outer diameter at the tube tip · I.D. = inner diameter at cannula tip · O-tube = outer cannula · I-tube = inner cannula · Length = length over the outer curve Length at mid radius = length of centre line from neck flange to tip