




FAHL® CRICOTHYROTOMY-SET

- ▶ To open up the airway in cases of emergency, i.e. in cases where no oral intubation is possible
- ▶ Set contains:
Cricothyrotomy cannula with trocar, flexible tube extension, tube holder, atraumatic suction catheter (Ch 10)
- ▶ The cricothyrotomy cannula with trocar consists of a straight cannula with ventilation connector (15 mm-connector/UNI) and an insertable trocar with integrated knife
- ▶ The plastic wings of the cannula limit the depth of penetration into the trachea and thereby reduce the risk of injuring the rear wall of the trachea.
- ▶ A „fixation bulge“ at the end of the cannula serves to avoid unintentional slipping of the cannula out of the trachea.
- ▶ If the trocar is removed, the tube extension, an anaesthesia bag or a ventilation connection (goose neck) can be connected via the straight 15 mm-connector (UNI) for ventilation.
- ▶ The tube holder enables the cannula to be fixed in place after the cricothyrotomy has been performed successfully.



STERILE EO Sterilized using ethylene oxide

⊗ Do not reuse

☀ Keep away from sunlight

☂ Keep dry

⊗ Do not use if package is damaged

📖 Consult instructions for use

CE 0482 CE marking identification number of the notified body

Product made in Germany

Order information:

FAHL® Cricothyrotomy-SET

REF 75001