



SPIRAFLEX® UNI / SPIRAFLEX® UNI CUFF

- ▶ Tracheostomy tube with adjustable neck flange to accommodate below or above average neck sizes, particularly for obese patients
- ▶ With cuff (REF 17503) or without cuff (REF 17504)
- ▶ Neck flange adjustable via reliable screw mechanism to avoid inadvertent sliding and dislocation of neck-flange
- ▶ A scale (5-mm increments) on the dorsal part of the tube allows exact and reproducible adjustment of the neck-flange
- ▶ With 15 mm connector (UNI) to connect to ventilation equipment or HMEs
- ▶ Set includes broad and adjustable tube holder, obturator and product ID card
- ▶ In addition to provide an atraumatic insertion of the tube the obturator allows suctioning of secretions. A hose connector with fingertip at the distal end of the obturator allows the connection of the suction line of a suction unit.



Order information:

SPIRAFLEX® UNI, SIZE 6-11



REF 17504

SPIRAFLEX® UNI CUFF, SIZE 6-11



REF 17503

When ordering, please state the size after the item number!

STERILE EO Sterilized using ethylene oxid

Do not resterilize

Keep away from sunlight

Single patient use

Consult instructions for use

Rx only CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner (in the US only).

CE 0482 CE marking identification number of the notified body

Product made in Germany

Size	O.D. neck flange mm	O.D. Tip mm	I.D. mm	Length mm	REF	REF with CUFF
6	8.4	8.4	6.0	adjustable to 52	17504-06	17503-06
7	10.2	10.2	7.0	adjustable to 68	17504-07	17503-07
8	11.2	11.2	8.0	adjustable to 88	17504-08	17503-08
9	12.2	12.2	9.0	adjustable to 117	17504-09	17503-09
10	13.2	13.2	10.0	adjustable to 118	17504-10	17503-10
11	14.2	14.2	11.0	adjustable to 123	17504-11	17503-11

O.D. Tip = Outer diameter of the outer tube at the tip of the cannula · O.D. neck flange = Outer diameter behind the neck flange · I.D. = Inner diameter at the tip of the cannula · Length = Length over the outer curve